



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mains Powered Vaccine Refrigerator 211 L

Product Code :
LBNY-0004-2100014



Description :

Vaccine storage capacity: 211L
Volume category: >= 120L
Refrigerant: R600a, 80 grams
Operating temperature: Min and Max rated ambient temperature: +5 °C and +43°C
UIFPF- Grade A
Prequalified Intermittent Mains Powered Refrigerator for storage of vaccines at +2 °C to +8°C.
Equipment operating voltage and frequency: 220-240 Volt and frequency 50Hz / 60Hz, 1ph.

Technical Specification :

Mains Powered Vaccine Refrigerator 211 L



Laboratory Instrument India