



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Mains Powered Combined Vaccine Refrigerator and Water  
pack Freezer 30L

**Product Code :**  
LBNY-0004-2100015



**Description :**

Volume category: Vaccine storage capacity: 30L

Refrigerant: R600a, (45 +40) grams

Operating temperature: Min and Max rated ambient temperature: +5 °C and +43°C

UIFPF- Grade A

PQS Prequalified Intermittent Mains Powered Combined Refrigerator and Water-Pack Freezer for storage of vaccines at +2 °C to +8°C and freezing of water packs.

220-240 Volt or 50Hz / 60Hz, 1ph

**Technical Specification :**

Mains Powered Combined Vaccine Refrigerator and Water pack Freezer 30L



**Laboratory Instrument India**