



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mains Powered Combined Vaccine Refrigerator and Water
pack Freezer 30L

Product Code :
LBNY-0004-2100015



Description :

Volume category: Vaccine storage capacity: 30L

Refrigerant: R600a, (45 +40) grams

Operating temperature: Min and Max rated ambient temperature: +5 °C and +43°C

UIFPF- Grade A

PQS Prequalified Intermittent Mains Powered Combined Refrigerator and Water-Pack Freezer for storage of vaccines at +2 °C to +8°C and freezing of water packs.

220-240 Volt or 50Hz / 60Hz, 1ph

Technical Specification :

Mains Powered Combined Vaccine Refrigerator and Water pack Freezer 30L



Laboratory Instrument India