



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Classroom DNA Electrophoresis Lab Station

**Product Code :**  
LBNY-0003-830008



**Description :**

Classroom DNA Electrophoresis Lab Station

**Technical Specification :**

Classroom DNA Electrophoresis Lab Station

Complete Electrophoresis system

Electrophoresis Apparatus-UV-transparent bottom, 2 comb positions, gel size 10 x 10 cm; made of UV and chemical resistant Plexiglas; Including safety plug with integrated power connections and comb with 12 pockets.

Electrophoresis Power supply-Stabilized power supply for both DNA and protein separation. Adjustable output voltage, timer with alarm function, two outputs for connecting of at least 2 electrophoresis chambers simultaneously.

Output voltage: 300 Volt (in 50 Volt stages),

Output current: 400 mA,

Power consumption: 60 Watt,

Input voltage: 100 – 240 Volt, AC.

Timer function: 1-999 min. with alarm function,

Gel visualization system

Gel electrophoresis kit- for 8 groups

---

# Laboratory Instrument India

