



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Mains Powered Combined Vaccine Refrigerator and Water  
pack Freezer 58 L

**Product Code :**  
LBNY-0004-2100016



**Description :**

Volume category: Vaccine storage capacity: 58 L

Energy consumption, stable running (kWh/24 hours): 1.91

Equipment operating voltage and frequency: 220-240 Volt 50Hz /60Hz, or 100-127 Volt 50Hz /60Hz, 1ph.

Operating temperature: Min and Max rated ambient temperature: +10 °C and +43°C

UIFPF- Grade A

Refrigerant: R600a, 68 grams

Prequalified Intermittent Mains Powered Combined Refrigerator and Water-Pack Freezer for storage of vaccines at +2 °C to +8°C and freezing of water packs.

**Technical Specification :**

Mains Powered Combined Vaccine Refrigerator and Water pack Freezer 58 L



**Laboratory Instrument India**