



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Mains Powered Vaccine Refrigerator 100 L

**Product Code :**  
LBNY-0004-2100017



**Description :**

Vaccine storage capacity: 100L  
Volume category:  $\geq 60$ , PQS Prequalified Intermittent Mains Powered Refrigerator for storage of vaccines at  $+2^{\circ}\text{C}$  to  $+8^{\circ}\text{C}$ .  
Operating temperature: Min and Max rated ambient temperature:  $+5^{\circ}\text{C}$  and  $+43^{\circ}\text{C}$   
UIFPF- Grade A  
Refrigerant: R600a, 80 grams  
Equipment operating voltage and frequency: 220-240 Volt and frequency 50Hz / 60Hz, 1ph.

**Technical Specification :**

Mains Powered Vaccine Refrigerator 100 L



**Laboratory Instrument India**