



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mains Powered Vaccine Refrigerator 100 L

Product Code :
LBNY-0004-2100017



Description :

Vaccine storage capacity: 100L
Volume category: ≥ 60 , PQS Prequalified Intermittent Mains Powered Refrigerator for storage of vaccines at $+2\text{ }^{\circ}\text{C}$ to $+8\text{ }^{\circ}\text{C}$.
Operating temperature: Min and Max rated ambient temperature: $+5\text{ }^{\circ}\text{C}$ and $+43\text{ }^{\circ}\text{C}$
UIFPF- Grade A
Refrigerant: R600a, 80 grams
Equipment operating voltage and frequency: 220-240 Volt and frequency 50Hz / 60Hz, 1ph.

Technical Specification :

Mains Powered Vaccine Refrigerator 100 L



Laboratory Instrument India