



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Mains Powered Vaccine Refrigerator 241 L

**Product Code :**  
LBNY-0004-2100018



**Description :**

Vaccine storage capacity: 241 L  
Volume category: >120L  
Energy consumption, stable running (kWh/24 hours): 2.21  
Operating temperature: Min and Max rated ambient temperature: +10°C and +43°C  
UIFPF- Grade A  
Refrigerant: R290, 145 grams  
Prequalified Intermittent Mains Powered Refrigerator for storage of vaccines at +2 °C to +8°C.  
Equipment operating voltage and frequency: 220-240 Volt, 50Hz/60Hz or 100-127 Volt 50Hz/60Hz, 1ph.

**Technical Specification :**

Mains Powered Vaccine Refrigerator 241 L



**Laboratory Instrument India**