



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mains Powered Vaccine Refrigerator 241 L

Product Code :
LBNY-0004-2100018



Description :

Vaccine storage capacity: 241 L
Volume category: >120L
Energy consumption, stable running (kWh/24 hours): 2.21
Operating temperature: Min and Max rated ambient temperature: +10°C and +43°C
UIFPF- Grade A
Refrigerant: R290, 145 grams
Prequalified Intermittent Mains Powered Refrigerator for storage of vaccines at +2 °C to +8°C.
Equipment operating voltage and frequency: 220-240 Volt, 50Hz/60Hz or 100-127 Volt 50Hz/60Hz, 1ph.

Technical Specification :

Mains Powered Vaccine Refrigerator 241 L



Laboratory Instrument India