



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mains Powered Vaccine Refrigerator 242 L

Product Code :
LBNY-0004-2100020



Description :

Vaccine storage capacity: 242 L

Volume category: >120L

Refrigerant: R600a, 130 grams

Energy consumption, stable running (kWh/24 hours): 0.638

Operating temperature: Min and Max rated ambient temperature: +5°C and +43°C

UIFPF- Grade A

PQS Prequalified Intermittent Mains Powered Refrigerator for storage of vaccines at +2 °C to+8°C.

Equipment operating voltage and frequency: 220-240 Volt 50Hz /60Hz, or 100-127 Volt 50Hz /60Hz, 1ph.

Technical Specification :

Mains Powered Vaccine Refrigerator 242 L



Laboratory Instrument India