



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Mains Powered Combined Vaccine Refrigerator and Water pack Freezer 120 L

**Product Code :**  
LBNY-0004-2100022



**Description :**

Volume category: 60L >60, Vaccine storage capacity: 120L  
Refrigerant: R600a, 57 grams refrigerator and 47 freezer compartments.  
Energy consumption, stable running (kWh/24 hours): 1.55  
Equipment operating voltage and frequency: 110-240 Volt 50Hz / 60Hz, 1ph.  
Operating temperature: Min and Max rated ambient temperature: +5°C and +43°C  
UIFPF- Grade A  
PQS Prequalified Intermittent Mains Powered Combined Refrigerator and Water-Pack Freezer for storage of vaccines at +2 °C to +8°C and freezing of water packs.

**Technical Specification :**

Mains Powered Combined Vaccine Refrigerator and Water pack Freezer 120 L



**Laboratory Instrument India**