



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mains Powered vaccine/water pack freezer 265 L

Product Code :
LBNY-0004-2100023



Description :

Gross volume-Freezer: 265L
Volume category: >150L
Refrigerant: R600a, 80 grams
Energy consumption, stable running (kWh/24 hours): 1.43
Equipment operating voltage and frequency: 220-240 Volt, and frequency 50Hz /60Hz, 1ph.
Operating temperature: Max rated ambient temperature: +43°C
PQS Prequalified Mains Powered Vaccine/Water-Pack Freezer for storage of vaccines at -15°C to -25°C and freezing of water packs.

Technical Specification :

Mains Powered vaccine/water pack freezer 265 L



Laboratory Instrument India