



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mains Powered Vaccine Refrigerator 145 L

Product Code :
LBNY-0004-2100028



Description :

Vaccine storage capacity: 145 L
Volume category: >120L
Energy consumption, stable running (kWh/24 hours): 0.69
Operating temperature: Min and Max rated ambient temperature: +5°C and +43°C
UIFPF- Grade A
R600a, 63 grams
Energy consumption, stable running (kWh/24 hours): 0.67
PQS Prequalified Intermittent Mains Powered Refrigerator for storage of vaccines at +2 °C to +8°C.
Equipment operating voltage and frequency.
220-240 Volt 50Hz /60Hz, or 100-127 Volt 50Hz /60Hz, 1ph.

Technical Specification :

Mains Powered Vaccine Refrigerator 145 L



Laboratory Instrument India

