



Email : sales@laboratoryinstrumentindia.com

Product Name :
Digital Scales

Product Code :
LBNY-0003-830010



Description :

Digital Scales

Technical Specification :

Digital scales
High precision sensor keep weighing stable and fast
Full range tare, Counting, Units choose.
RS 232 Interface can connect the computer and printer.
Rechargeable battery.
Dual display and dual line display can be choose.
Automatic calibration
Range: Up to 200 g
Resolution: $\pm 0,01$ g
Calibration, TARA, units conversion,
Power: 220V/50Hz,



Laboratory Instrument India

