



Email : sales@laboratoryinstrumentindia.com

Product Name :
Solar Direct Drive Vaccine Refrigerator with no battery 220

Product Code :
LBNY-0005-2100017

L



Description :

Vaccine storage capacity: 220 L

Volume category: >120L

Solar radiation reference period of 3.5 kWh/m²/day and hot zone (+43°C) test temperatures.

Refrigerant: R600a, 65 grams

PQS prequalified Solar direct drive refrigerator with no battery for storage of vaccines at +2 °C to +8°C in areas without any electricity, or where there is less than eight hours of reliable electricity over a typical day.

Operating temperature: Min and Max rated ambient temperature: +5 °C and +43°C

UIFPF: Grade A

Technical Specification :

Solar Direct Drive Vaccine Refrigerator with no battery 220 L



Laboratory Instrument India