



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Digital Multipurpose scale

**Product Code :**  
LBNY-0003-830011



**Description :**

Digital Multipurpose scale

**Technical Specification :**

Digital scales

It can be used in multiple places as per your choice and about how you want to use it. This is the best weighing machine in the market.

Multipurpose scale. Used for weighing all general items like grocery, dry fruits, industrial products etc. Can also be used for weighing costly medicines and chemicals also.

With support for running on direct AC using a standard adaptor . This advanced version not only can measure but also has a feature to count in units.

Auto Shut Off Counting Function: Low Battery Indication this equipped with a high precision strain Gauge Sensor System. Has TARE option to weigh different ingredients in the same container.

Automatic calibration

Range: Up to 500 g

Resolution:  $\pm 0,1g$

Calibration, TARA, units conversion,

Power: 220V/50Hz,

---

# Laboratory Instrument India

