



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Health Center Kit

**Product Code :**  
LBNY-0005-2100021



**Description :**

An integrated kit linking SDD refrigerators to a solar generator, which uses excess electricity from solar panels to power lighting, mobile phones, ventilation and other applications.

- 1 x mobile ceiling fan,
- 2 x USB outputs (5W each),
- 1 x 12V/20W socket outlet,
- 3 x keys for the locker,
- 1 x rechargeable 27Ah battery with a life cycle of 5 years,
- 2 x rechargeable LED lights with 3 settings and an autonomy of 35 hours.

**Technical Specification :**

Health Center Kit



**Laboratory Instrument India**