



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Solar Direct Drive Combined Vaccine Refrigerator and  
Waterpacks Freezer with no battery 120 L

**Product Code :**  
LBNY-0005-2100025



**Description :**

Vaccine storage capacity: 120 L

Volume category:  $\geq 120L$

Autonomy time: Appliance's autonomy at +43 °C ambient: 82 hours 5 minutes

Refrigerant: R600a, 57 grams / Freezer: 47 grams

Operating temperature: +5 °C to +43°C

UIFPF: Grade A

PQS prequalified Solar direct drive combined refrigerator and freezer with no battery for storage of vaccines at +2 °C to +8°C and freezing of ice packs

Solar radiation reference period of 3.5 kWh/m<sup>2</sup>/day and hot zone (+43°C) test temperatures.

**Technical Specification :**

Solar Direct Drive Combined Vaccine Refrigerator and Waterpacks Freezer with no battery 120 L



**Laboratory Instrument India**

