



Email : sales@laboratoryinstrumentindia.com

Product Name :
Solar Direct Drive Combined Vaccine Refrigerator and
Waterpacks Freezer with no battery 120 L

Product Code :
LBNY-0005-2100025



Description :

Vaccine storage capacity: 120 L

Volume category: $\geq 120L$

Autonomy time: Appliance's autonomy at $+43\text{ }^{\circ}C$ ambient: 82 hours 5 minutes

Refrigerant: R600a, 57 grams / Freezer: 47 grams

Operating temperature: $+5\text{ }^{\circ}C$ to $+43\text{ }^{\circ}C$

UIFPF: Grade A

PQS prequalified Solar direct drive combined refrigerator and freezer with no battery for storage of vaccines at $+2\text{ }^{\circ}C$ to $+8\text{ }^{\circ}C$ and freezing of ice packs

Solar radiation reference period of $3.5\text{ kWh/m}^2/\text{day}$ and hot zone ($+43\text{ }^{\circ}C$) test temperatures.

Technical Specification :

Solar Direct Drive Combined Vaccine Refrigerator and Waterpacks Freezer with no battery 120 L



Laboratory Instrument India

