



Email : sales@laboratoryinstrumentindia.com

Product Name :
Digital Analytical Scales

Product Code :
LBNY-0003-830012



Description :

Digital Analytical Scales

Technical Specification :

Digital Analytical scales
High precision sensor keep weighing stable and fast
Full range tare, Counting, Units choose.
RS 232 Interface can connect the computer and printer.
Rechargeable battery.
Dual display and dual line display can be choose.
Automatic calibration
Range: Up to 200 g
Resolution: $\pm 0,001$ g
Calibration, TARA, units conversion,
Power: 220V/50Hz,



Laboratory Instrument India

