



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Digital Analytical Scales

**Product Code :**  
LBNY-0003-830012



**Description :**

Digital Analytical Scales

**Technical Specification :**

Digital Analytical scales  
High precision sensor keep weighing stable and fast  
Full range tare, Counting, Units choose.  
RS 232 Interface can connect the computer and printer.  
Rechargeable battery.  
Dual display and dual line display can be choose.  
Automatic calibration  
Range: Up to 200 g  
Resolution:  $\pm 0,001g$   
Calibration, TARA, units conversion,  
Power: 220V/50Hz,



**Laboratory Instrument India**

