



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
30DTR,Int. Sensor,Refr

**Product Code :**  
LBNY-0006-210001



**Description :**

The device shows the actual temperature, all alarm violations over the previous 30 days (on a rolling basis), the daily minimum and maximum temperature of the last 30 days, and the time duration of any violation. with internal sensor, is an electronic 30-Day Temperature Recorder (30DTR), with factory-programmed alarms and visual display for monitoring storage conditions in vaccine refrigerators (+2°C to +8°C) over a 30-day period.

Further data analytics such as daily min/max and average temperatures, time out of range, alarm history etc. can be generated via the free HETL-01 Data Management Software.

USB interface and can be connected to a PC or Mac without any installation software.

When connected, will automatically generate a 30-day raw data report in CSV file format.

**Technical Specification :**

30DTR,Int. Sensor,Refr



**Laboratory Instrument India**