




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
30DTR Accessory USB Cradle

**Product Code :**  
LBNY-0006-2100014



#### **Description :**

The USB Interface Cradle is supplied with a 1-meter micro-USB to USB-A cable for connection to a PC and can be used to download data from multiple 30DTRs consecutively. The USB Interface Cradle serves as the physical data connection between a 30-Day Temperature Recorder (30DTR) and a PC for temperature data download.

#### **Technical Specification :**

##### **Hardware Specifications:**

Operating temperature: -20°C to +60°C

Material: Polycarbonate, ABS

User interface: Red and green LEDs

USB standard: 2.0

Power Source: Powered from PC

IP rating: IP50

Connection: Via micro-USB to USB-A cable (length 1 meter)

Unit Dimensions (LxWxD): 8.5 x 7.1 x 4.5cm

Unit Weight: 74g.



**Laboratory Instrument India**

