



Email : sales@laboratoryinstrumentindia.com

Product Name :
IFI Tag

Product Code :
LBNY-0006-2100015



Description :

The device indicates via a visual alarm when an exposure to temperatures below -0.5°C for longer than 60 minutes has occurred.

An electronic Irreversible Freeze Indicator (IFI) with an internal sensor, factory-programmed alarm, and visual display for indicating a freeze event during transport or storage of freeze-sensitive vaccines.

The device needs to be activated prior to use and once an alarm has been triggered, cannot be re-used.

Technical Specification :

Hardware Specifications:

Operating temperature: -25°C to 60°C

Alarm: Visual, display showing either a tick for "OK" or a cross for "Alarm"

User interface: LCD, with start button switch

Sensor: Internal, non-replaceable sensor

IP rating: IP64

Memory: Not applicable

Data download: Not available

Mounting: Mounting lug

Label colour: Blue

Resolution: 0.1°C

Accuracy: $\pm 0.5^{\circ}\text{C}$ within the range -10°C to $+25^{\circ}\text{C}$

Battery type: Non-replaceable lithium coin cell

Shelf life: 42 months from date of manufacture (includes activated life)

Activated life: 36 months

Unit Dimensions (LxWxD): 7.15 x 3.3 x 0.86cm
Unit weight: 19g.



Laboratory Instrument India