



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Tool kit, f. ICT technicians, 110 VAC eqpt

**Product Code :**  
LBNY-0008-210001



**Description :**

Tool kit composed of essential hand-tools for installation and maintenance of telecommunications systems and solutions (radio, VSAT, data-communications, etc) at country offices.

**Technical Specification :**

Tool kit composed of essential hand-tools (screwdrivers, wrenches, pliers, spanners, coaxial cable crimping tools, two soldering irons (electrical and butane-gas operated)), drill bits, electrical and battery drills, measuring instruments (SWR-meter, VOM Multimeter), flashlights, compass and inclinometer (for VSAT alignment), measuring tape, knives, cable ties, hacksaw, Allen keys, antistatic wrist-band and various other useful tools



**Laboratory Instrument India**