



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tool kit, f. ICT technicians, 220 VAC eqpt

Product Code :

LBNY-0008-210002



Description :

Tool kit composed of essential hand-tools for installation and maintenance of telecommunications systems and solutions (radio, VSAT, data-communications, etc.) at country offices.

Technical Specification :

Tool kit composed of essential hand-tools (pliers, spanners, screwdrivers, wrenches, coaxial cable crimping tools, two soldering irons (electrical and butane-gas operated)), electrical and battery drills, drill bits, measuring instruments (SWR-meter, VOM Multimeter), flashlights, compass and inclinometer (for VSAT alignment), hacksaw, Allen keys, measuring tape, knives, cable ties, antistatic wrist-band and various other useful tools.



Laboratory Instrument India