



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
HF Radio Base Kit

**Product Code :**  
LBNY-0009-210001



**Description :**

HF base station radio kit; PSU, power cables w/fuse, radio w/handset, RF unit, desk console, base antenna, coaxial cable, lightning arrestor  
1.6-30MHz, 125W PEP, 400 channels, USB/LSB, AM, AFSK, scanning  
Data ready with additional cooling fans for interconnect with data modem  
Standard 4/6-digit selective calling, emergency calling, GPS enabled, ALE/CALM, optional telephone interconnect.

**Technical Specification :**

HF Radio Base Kit



**Laboratory Instrument India**