



Email : sales@laboratoryinstrumentindia.com

Product Name :
HF Radio Base Kit

Product Code :
LBNY-0009-210001



Description :

HF base station radio kit; PSU, power cables w/fuse, radio w/handset, RF unit, desk console, base antenna, coaxial cable, lightning arrestor
1.6-30MHz, 125W PEP, 400 channels, USB/LSB, AM, AFSK, scanning
Data ready with additional cooling fans for interconnect with data modem
Standard 4/6-digit selective calling, emergency calling, GPS enabled, ALE/CALM, optional telephone interconnect.

Technical Specification :

HF Radio Base Kit



Laboratory Instrument India