



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
HF Radio Mobile Kit

**Product Code :**  
LBNY-0009-210002



**Description :**

HF mobile station radio kit; radio w/handset, RF unit, junction box, power cables w/fuse, mobile antenna w/two whips, coaxial cable, control cable.

1.6-30MHz, 125W PEP, 400 channels, USB/LSB, AM, AFSK, scanning.

Data ready with additional cooling fans for interconnect with data modem

Standard 4/6-digit selective calling, emergency calling, GPS enabled, ALE/CALM, optional telephone interconnect.

**Technical Specification :**

HF Radio Mobile Kit



**Laboratory Instrument India**