




Email : sales@laboratoryinstrumentindia.com

Product Name :
HF Data Modem

Product Code :
LBNY-0009-210003



Description :

Data mode: Selective repeat ARQ Protocol, 2400 bps-16 channel QPSK modulation.
Transceiver Interface: 9600 baud RS232; TR.29 based AT port.
Primary Power: 13.5V DC nominal (250 mA max), 9-16V DC operating range.
Link establishment mode: Proprietary link establishment 80 baud CHIRP.
Temperature: 0°C to 55°C.
Data compression: In-built data compression.
Rate of data transfer: Up to 6000 pbs (compressed), Up to 1475 bps (uncompressed).
Size: 210mm W x 240mm D x 65mm H (including rear connectors).
Weight: 1,8kg.

Technical Specification :

HF Data Modem



Laboratory Instrument India