




Email : sales@laboratoryinstrumentindia.com

Product Name :
Portable Antenna Mast, Aluminum, 12m

Product Code :
LBNY-0009-210006



Description :

Suitable for preliminary emergency response set-up or small size field office installations. It is not recommended as a permanent solution or for large and heavy antennas. Complete portable antenna mast supporting light-weight HF and VHF antennas for temporary installations.

Technical Specification :

Height: 12m pulled out
Guy radius: 6m
Portable antenna mast kit: four-section collapsible mast, rope guys, erection kit, ground plate, masthead gibbet and halyard pulley
Wind loading: 12kg overhead load, 144km/h (90mph)
Weight: 28.4kg.



Laboratory Instrument India