



Email : sales@laboratoryinstrumentindia.com

Product Name :
Register, A4, squared, 80 pgs

Product Code :
LBNY-0010-210001



Description :

Register, A4, squared, 80 pages, pack of 10

Technical Specification :

Finishing: Saddle stitch (stainless steel) along the 297 mm and trim to size.

Printing: One color Cyan blue and water based glossy varnish recto only cover and down the free side of the back cover, in white.

Register, A4 size (210 x 297mm), white wood-free paper of 70gsm suitable for both pen and pencil, 40 sheets/80 pages. Pages printed with 5mm blue grey squares with no margin.

Cover: Glossy finish, 270gsm offset paper.



Laboratory Instrument India