



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Register, A4, squared, 80 pgs

**Product Code :**  
LBNY-0010-210001



**Description :**

Register, A4, squared, 80 pages, pack of 10

**Technical Specification :**

Finishing: Saddle stitch (stainless steel) along the 297 mm and trim to size.

Printing: One color Cyan blue and water based glossy varnish recto only cover and down the free side of the back cover, in white.

Register, A4 size (210 x 297mm), white wood-free paper of 70gsm suitable for both pen and pencil, 40 sheets/80 pages. Pages printed with 5mm blue grey squares with no margin.

Cover: Glossy finish, 270gsm offset paper.



**Laboratory Instrument India**