



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Paper, plain, 70cm wide, 25m ROLL

**Product Code :**

LBNY-0010-210005



**Description :**

The composition of the paper should be homogenous with no inclusions/alien particles and be suitable for a range of art activities (painting, colouring, sticking etc.)

Roll of white, plain (un-ruled) paper, 70 cm wide and 25 m long, preferably with an internal cardboard tube strengthener. Paper to be 75% recovered paper fibres (post and pre-consumer waste) and of 80 gsm.

**Technical Specification :**

Paper, plain, 70cm wide, 25m ROLL



**Laboratory Instrument India**