



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Micro Centrifuge

**Product Code :**  
LBNY-0003-830018



**Description :**

Micro Centrifuge

**Technical Specification :**

MicroCentrifuge

For quickly spinning down and mixing solutions during DNA electrophoresis and PCR experiments

Speed 6,000 rpm/2,000 x g

Safe on/off switch

Starts and stops in seconds

Capacity for 6 x 1.5/2.0 ml tubes



**Laboratory Instrument India**