



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Paint, Chalkboard, Black

**Product Code :**  
LBNY-0010-2100012



**Description :**

Paint, blackboard, black, non-inflammable, non-toxic, water-based, tin of 500ml, with spring locked lid to prevent spillage during transport.

The labelling on the tin provides clear instructions for use; instructions on how to clean the brush after painting, as well as storage conditions (frost-free).

It is essential that the opening of the tin is 65mm diameter, enough to accommodate the Paint Brush (width 50-60mm).

**Technical Specification :**

Paint, Chalkboard, Black



**Laboratory Instrument India**