



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Bag, Biodegradable PAC-20

**Product Code :**  
LBNY-0010-2100017



www.laboratoryinstrumentindia.com

**Description :**

Biodegradable carrier bag, printed with 4 colors, pack of 20.

**Technical Specification :**

Biodegradable carrier bag with handle.

Material: 100% biodegradable.

Color: white

GSM:50GSM

Print: Up to 4 process colors. 70% coverage on both sides

The "zip lock" seal on the bag is secured and strong enough to carry up to 1000g of school materials. Books of A4 size are easily put into the bag, sealed, and removed.

Dimensions: 430mm length x 290mm width. The bag is heat sealed on all sides including the handle which is double layered for added strength. active length (excluding the sealed edges and handle) (length below the base of the seal) of 350 mm, length above the top of the seal of 80 mm, width of 282 mm and thickness of 0.15 mm (150 microns).The bag must contain a "zip lock" seal below the handle.



**Laboratory Instrument India**

