



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Bag, Plastic, Interlock Seal/PAC-20

**Product Code :**

LBNY-0010-2100018



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

100% recycled LDPE Plastic carrier bag, self sealing with interlock seal type closure on 1 short side, securely sealed along the 2 long sides with a cleanly punched handle.

Material: 100% recycled LDPE

Color: white

Dimensions: active length (length below the base of the seal) of 350mm, length above the top of the seal of 80mm, width of 282mm and thickness of 0.15mm (150 microns).

**Technical Specification :**

Bag, Plastic, Interlock Seal/PAC-20



**Laboratory Instrument India**