



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bag, Plastic, Interlock Seal/PAC-20

Product Code :
LBNY-0010-2100018



www.laboratoryinstrumentindia.com

Description :

100% recycled LDPE Plastic carrier bag, self sealing with interlock seal type closure on 1 short side, securely sealed along the 2 long sides with a cleanly punched handle.

Material: 100% recycled LDPE

Color: white

Dimensions: active length (length below the base of the seal) of 350mm, length above the top of the seal of 80mm, width of 282mm and thickness of 0.15mm (150 microns).

Technical Specification :

Bag, Plastic, Interlock Seal/PAC-20



Laboratory Instrument India