



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Box, metal, lockable, for storage

**Product Code :**

LBNY-0010-2100019



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

**Metal Box:**

Made of high density, durable lightweight aluminum. Robust, lockable, stackable. Water proof, dustproof and impact resistant.

Withstand incredible abuse and extremely harsh environments without chipping, denting, or losing their finish.

Code Padlocks, set of 2, 2 sets are provided with the box:

The padlock consists of shackle, a dial, locking parts, and means for combination changing and combination scrambling when locked. Locks are fixed inside the cover of the box.

Material, recycled aluminum, thickness 1mm.

Internal dimensions: 764x574x220mm

External dimensions: 790x590x245mm

2 Strap for cover nylon

2 Locking front side lock for code padlocks.

**Technical Specification :**

Box, metal, lockable, for storage



**Laboratory Instrument India**

