



Email : sales@laboratoryinstrumentindia.com

Product Name :

Decal, equip., small, sheet of 10

Product Code :

LBNY-0012-210001



www.laboratoryinstrumentindia.com

Description :

Decal, equipment, small, 65x25mm (sheet of 10)

Material: White gloss PVC, with permanent adhesive; 0.1mm thickness; lifespan 2-3 years from date of delivery; Kraft liner 130 g.

Color: Background color; cyan blue. The logo printed in white, positioned in the middle of the decal, running 80% of the length.

Size: Die cut in size 65 x 25 mm, sheet of 10 decals (2 x 5 pieces).

Decal, permanently adhesive, size: 65x25mm, sheet of 10 decals (2 x 5 pieces).

Technical Specification :

Decal, equip., small, sheet of 10



Laboratory Instrument India