



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Water Distiller

**Product Code :**  
LBNY-0003-830020



**Description :**

Water Distiller

**Technical Specification :**

Water distiller

Capacity; 5 l/h

Voltage: 220V

Stainless steel with carboy with spigot and screw closure

Made of strong, scratch-resistant polypropylene with excellent chemical resistance and capacity of min 10 L



**Laboratory Instrument India**