



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Decal, vehicle, small, 693x196mm

**Product Code :**

LBNY-0012-210002



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Material: White gloss PVC, with permanent adhesive; 0.1mm thickness; lifespan 2-3 years from date of delivery; Kraft liner 130 grams

Colour: Background colour: cyan blue. The logo printed in white, positioned in the middle of the decal, running 80% of the length.

Decal, vehicle, permanently adhesive, 693 x 196 mm

Size: Rectangular, length 693 x 196 mm, with corners in a 90-degree angle.

**Technical Specification :**

Decal, vehicle, small, 693x196mm



**Laboratory Instrument India**