



Email : sales@laboratoryinstrumentindia.com

Product Name :

Decal, rect, 205x105mm/RL-50

Product Code :

LBNY-0012-210003



www.laboratoryinstrumentindia.com

Description :

Decal, equipment, permanently adhesive

Roll of 50

Material: White gloss PVC with permanent adhesive, 0.1mm thickness, lifespan 2-3 years from date of delivery, Kraft liner 67 grams.

Base color: 100% Cyan blue. The logo are printed in white, and positioned in the middle of the decal, running 80% of the length.

Decal, rectangular, 205 x 105mm, roll of 50 decals

Size: Rectangular, length 205 mm x height 105 mm, with corners in a 90-degree angle.

Technical Specification :

Decal, rect, 205x105mm/RL-50



Laboratory Instrument India