



Email : sales@laboratoryinstrumentindia.com

Product Name :
Decal, Round, Diameter 205mm

Product Code :
LBNY-0012-210004



www.laboratoryinstrumentindia.com

Description :

Decal, round, permanently adhesive, on square sheet.

Material: White gloss PVC, with permanent adhesive; 0.1mm thickness; lifespan 2-3 years from date of delivery; Kraft liner 130 grams.

Color: Background colour; cyan blue. The printed in white, positioned in the middle of the decal, running 80% of the diameter. Outside format: White background, with instructions printed in cyan blue lettering.

Size: Die-cut to diameter 205 mm, delivered in squares with instructions printed in cyan blue lettering on the white outside format.

Technical Specification :

Decal, Round, Diameter 205mm



Laboratory Instrument India