



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Decal, Vehicle, Large, 1000x280mm

**Product Code :**

LBNY-0012-210005



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Decal, vehicle, permanently adhesive, 1000 x 280mm

Size: Rectangular, length 1000 mm x height 280 mm, with corners in a 90-degree angle.

Material: White gloss PVC, with permanent adhesive; 0.1mm thickness; lifespan 2-3 years from date of delivery; Kraft liner 130 grams.

Color: Background color: cyan blue. The logo printed in white, positioned in the middle of the decal, running 80% of the length.

**Technical Specification :**

Decal, Vehicle, Large, 1000x280mm



**Laboratory Instrument India**