



Email : sales@laboratoryinstrumentindia.com

Product Name :
Insignia, Cloth, 85 x 41mm

Product Code :
LBNY-0012-210006



www.laboratoryinstrumentindia.com

Description :

Material: Cotton/polyester blends.
The edges of the insignia overlock-stitched or heat cut, to avoid fraying.
Color: Cyan blue, with a woven horizontal logo in white, centered vertically.
Insignia, cotton/polyester blend, rectangular, 85x41mm.

Technical Specification :

Insignia, Cloth, 85 x 41mm



Laboratory Instrument India