



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Insignia, Cloth, 85 x 41mm

**Product Code :**  
LBNY-0012-210006



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Material: Cotton/polyester blends.  
The edges of the insignia overlock-stitched or heat cut, to avoid fraying.  
Color: Cyan blue, with a woven horizontal logo in white, centered vertically.  
Insignia, cotton/polyester blend, rectangular, 85x41mm.

**Technical Specification :**

Insignia, Cloth, 85 x 41mm



**Laboratory Instrument India**