



Email : sales@laboratoryinstrumentindia.com

Product Name :
Round Metal Badge, 55 mm Diameter

Product Code :
LBNY-0012-210007



www.laboratoryinstrumentindia.com

Description :

Material: Metal base
Badge, metal, round, cyan blue
Size: Diameter 55mm
Position of Logo: Logo vertically centered, running 80% of the diameter of the badge.
Color: Cyan blue.
Logo height: 7mm.
Logo length: 47mm.

Technical Specification :

Round Metal Badge, 55 mm Diameter



Laboratory Instrument India