



Email : sales@laboratoryinstrumentindia.com

Product Name :
Armband, Cyan Blue Cloth

Product Code :
LBNY-0012-210008



www.laboratoryinstrumentindia.com

Description :

Armband, cloth.
Material - 220 gsm 100% washed cotton twill.
The material is color-fast, durable and easy to wash.
Adjustable Velcro closure consisting of two horizontal Velcro strips;
length 120mm x width 30mm; sewn along the outer side of the armband.
A third vertical Velcro strip sewn on the underside of the armband,
height 80mm x width 30mm
Color: 100% Cyan Blue
The colour is durable and does not bleach easily in strong sunlight.
There is no leaching of colour when the armband comes into contact with water.
Size: length 420mm x height 80mm.

Technical Specification :

Armband, Cyan Blue Cloth



Laboratory Instrument India

